


## IEXdata structure



## IEXdata → data class to handle loading of mda and EA files

```
Mydata=IEXdata(*scanNum, dtype='mda', **kwargs)
```

**\*scanNum**

- scanNum → single scan
- first, last → every scan from first to last inclusive
- first, last, countby → every nth scan from first to last inclusive
- last = inf → to load to the end of the directory
- [scan1, scan2, scan3] → series of scan

**dtype** → mda by default (EA, EA\_nc, nData)

**\*\*kwargs**

- path
- prefix
- suffix
- nzeros
- overwrite True/False; True → reloads the data, False → skips already loaded data

`mydata.update(*scanNum, **kwargs)` → loads additional scans to the dictionary (uses previous dtype and file info (path, prefix, suffix, nzeros))

`mydata.info()` → prints the loaded scans and current variable values

`remove(mydata.mda, * scanNum)` → remove scan from dictionary

`mydata.save(filename, filepath)` → to save a data set as hdf5

`reload_IEXdata(filename, filepath)` → to reload saved data

### nmda → data class for mda files

```
mydata= (scanNum*)
```

- `mydata.mda` → dictionary of mda scans
- `mydata.mda[scanNum]` → dictionary of all detector data from scanNum
- `myddata.mda[scanNum].det[detNum]` → pynData object
- `mydata.mda[scanNum].header` → IEX header object

### nEA → data class for EA files

```
mydata=IEXdata(scanNum*, dtype='EA')
```

- `mydata.EA` → dictionary of EA scans
- `mydata.EA[scanNum]` → pynData object with addition  
pynData\_ARPES vars
- `mydata.EA[scanNum].header` → IEX header object

# Loading data

## mda scans

Loading for the first time:

```
mydata = IEXdata(scanNum) → single scan  
mydata = IEXdata(first, last) → every scan from first to last inclusive  
mydata = IEXdata(first, last, countby) → every nth scan from first to last inclusive  
mydata = IEXdata([scan1, scan2, scan3]) → series of scan
```

**Note:** last = inf → to load to the end of the directory

Adding scans (uses the same scanNum syntax as above):

```
mydata.update(first, inf, overwrite=False) → loads all unloaded scans (overwrite=True is the default and will reload the already loaded data)
```

## EA scans

Loading for the first time (uses the same scanNum syntax as above):


```
mydata = IEXdata(scanNum, dtype='EA') → for new EA .h5 format  
mydata = IEXdata(scanNum, dtype='EA_nc') → for old EA .nc format
```

Adding scans (uses the same scanNum commands):

```
mydata.update(first, inf, overwrite=False) → loads all unloaded scans (overwrite=True is the default and will reload the already loaded data)
```

**Note:** update uses that last dtype for load unless otherwise specified

## nmda structure



# nmda

`nmda.fpath` → `pnData_ARPES` object

`nmda.header` → header object

`nmda.detAll` → print list of detectors `detNum: (PV, description, units)`

`nmda.det[detNum]` → `nData` object of `detNum` data (`.data` => data array)

`nmda.det[detNum].pv` → detector: (PV, description, units)

`nmda.posAll` → print list of positioners `posNum: (PV, description, units)`

`nmda.posx[posNum]` → `nData` object of `posNum` readback values (`.data` => data array)

`nmda.posx[posNum].pv` → positioner: (PV, description, units)


`nmda.posy[posNum]`, `nmda.posz[posNum]`, `nmda.post[posNum]` → for higher dimensional data

`nmda.setscalePos('x', PosNum)` → sets the 'x' / 'y' / 'z' scale for all detectors to the Positioner with index `Pos`

`nmda.setscalePos('x', DetNum)` → sets the 'x' / 'y' / 'z' scale for all detectors to scaling from values of `DetNum`

`nmda.setscaleIndex('x')` → sets the 'x' / 'y' / 'z' scale for all detectors index/point num

## nEA structure



## nEA / ndata\_ARPES

### nEA

nEA → pnData\_ARPES object

nEA.fpath → full path to the original file

nEA.mdafname → mda file name

nEA.header → IEX header object

### Ndata\_ARPES

nEA.EDC → pnData object of angle integrated data

nEA.MDC → pnData object of energy integrated data

nEA.hv → photon energy (var or array)

nEA.wk → analyzer work function (var or array)

nEA.thetaX → polar angle (var or array)

nEA.thetaY → other angle (var or array)


nEA.KEScale → array of original KE scaling

nEA.angScale → array of original detector angle scaling

nEA.angOffset → offset in original detector angle

nEA.slitDir → analyzer slit direction

## pynData structure





# pynData

`nData.data` → data array

`nData.scale[ 'x' ]` → data array for 'x' / 'y' / 'z' scale

`nData.unit[ 'x' ]` → units for 'x' / 'y' / 'z' scale

`nData.extras` → dictionary for meta data

`nData.info()` → prints data array shape and axis info

`nData.updateAx( 'x', NewscaleArray, 'Scale_units' )` → function for changing 'x' / 'y' / 'z' scale

`nData.shiftScale( 'x', oldValue, newValue )` → function for changing 'x' / 'y' / 'z' scale by calculating an offset from oldValue and newValue

`nData.save()` → saves pynData object as an hdf5 file `load_nData(fname, fdir='')` **Note: only saves nData.data, nData.scale and nData.extras**


`load_nData(fname, fdir='')` → load nData saved via .save()

# Simple plotting examples using pynData


## Niceplot(pynData\_object) simple 1D and 2D plotting

```
[30]: niceplot(data.mda[221].det[36],label="H")
      niceplot(data.mda[223].det[36],label="V")
      plt.legend()
```

```
[30]: <matplotlib.legend.Legend at 0x7f36d269feb8>
```



```
[70]: 1 niceplot(data.mda[31].det[15])
```




## plot2D(pynData\_object) plotting with linecuts

```
[79]: 1 plot2D(data.mda[31].det[15],xWidthPix=2)
```

```
10
```


```
[79]: (<Macros_29id.pynData.pynData.nData at 0x7f376482e850>,  
<Macros_29id.pynData.pynData.nData at 0x7f3756468790>)
```



## plot3D(pynData\_object)


```
[190]: plot3D(FM,xCen=496.75,xWidthPix=10)
```

```
[190]: (<pynData.nData at 0x7f52264138d0>,  
<pynData.nData at 0x7f52264139e8>,  
<pynData.nData at 0x7f5226413a20>,  
<pynData.nData at 0x7f52265a7b70>)
```



IEX specific plotting


# Interactive plotting



# Python plotting

```
[49]: #plt.plot(x,y)
plt.plot(data.mda[221].det[36].scale['x'],data.mda[221].det[36].data,label="H")
plt.plot(data.mda[223].det[36].scale['x'],data.mda[223].det[36].data*1.125+0.125,label="V")
plt.legend()
```

[49]: <matplotlib.legend.Legend at 0x7f36d1fd7550>



## Saving experiment: IEX\_nData hdf5 structure

h5 file

- group: mda => contains all the mda scans
  - group: mda\_scanNum
    - .attrs -> fpath (original file path)
  - group: header
    - .attrs for each header attribute
  - group: det
    - group: det\_1 .... -> nData with subgroups for each detector
  - group: posx, posy ...
    - group: posx\_1 .... -> nData with subgroups for each positioner
- group: EA => contains all the EA scans
  - group: EA\_scanNum -> nData with subgroups
    - group: EDC -> nData with subgroups
    - group: MDC -> nData with subgroups
    - .attrs: hv, wk, thetaX, thetaY, angOffset, slitDir
    - dataset: KEscale, angScale
    - .attrs fpath
    - .attrs mdafname
  - group: header
    - .attrs header info

- group: nDataObject
  - subgroup: data
    - dataset -> data array
  - subgroup: scale
    - dataset -> scale array
  - subgroup: unit
    - .attrs -> unit

Common command examples

```
f= h5py.File("/Users/jmcchesn/Documents/NoteBooks/Data/IEXmydata2.h5", 'r')
```

```
for k in f.keys():  
    print(k)
```

```
for k in f.attrs.keys():  
    print('{} => {}'.format(k, f.attrs[k]))
```

**View the data with python**





## **ARPES data files**

images: EA\_0001.h5

motor scans: ARPES\_0001.mda

